

Exp. 1 (Text #3): The Molecular Sieve Zeolite X

Lab Work 1/24

Report Due 2/6



What is a zeolite?

- A microporous solid, containing pores or channels in its structure that can accommodate guest molecules
- An aluminosilicate
 - Framework stoichiometry: $(\text{Si}, \text{Al})_n \text{O}_{2n}$
 - Si or Al atoms are tetrahedrally coordinated to bridging O's ("vertex-sharing" tetrahedra)
 - Cations (e.g., Na^+ , K^+) *required for charge balance*




Zeolite Structure

Several aluminosilicate structures are based on a **truncated octahedron** with stoichiometry $M_{24}O_{36}$ (where $M = Si, Al$), also called the **sodalite** or β cage:



Zeolite A
(showing connectivity of "octahedra")




Zeolite X

Zeolite Structure, continued

Cations occupy numerous sites within the framework, and help to determine the size of the pores (α - or supercage).

– Also influenced by Si/Al ratio



We will use Na^+ to balance charge, so the **hydrated sodium ion** helps determine pore size.

How would pore size change if K^+ were used, instead?

K^+ is larger; pore size would be smaller.

Applications of Zeolites

- Molecular sieves (separation by size):
 - Desiccants/Adsorbents
- Ion exchange
 - Water softening
- Catalysis
 - Introduction of transition-metal ions affords numerous sites for catalytic reactions

NaX Synthesis and Ion-Exchange:

Synthesis of NaX:



$n = 9.6 - 12$ for X-type zeolites; For us, $n = \mathbf{10.7}$ (pore size = 7.4 Å)

Completed by mixing prepared solutions of sodium aluminate and sodium silicate (Solutions 1 and 2 in text)

Characterization by IR spectroscopy (1/25): See Balkus, K. J. et al. *J. Chem. Educ.* **1991**, 68, 875-877 for published spectra.

Ion-Exchange Reaction (1/25):



What is the specific ion-exchange process that occurs here?

Uptake of 1 Co²⁺ results in release of 2 Na⁺ ions


NaX Synthesis: Procedural Notes and Tips

- You will work in pairs on this experiment.
- We will perform the experiment at **50% scale**.
- Next week, we will complete Part A to the stopping point mentioned in the text (filtering NaX crystals and leaving them in your drawer to dry).
- Our aim is to allow 2 hours for reaction in the oven, ideally starting around 3:00. We cannot begin heating until everyone is ready.
 - Make water bath immediately and start heating (watch temp as directed)
 - Work on Solutions 1 and 2 simultaneously
 - Note that the specified masses are not very precise (e.g., 1.2 g). **Don't** waste time trying for 1.200 g; just record the exact mass you obtain.
- After cooling, you will suction-filter your product using a Buchner funnel and filter paper. Wash the crystals with ~3 portions of water and continue suction as long as possible.
- Be careful when removing crystals from filter paper.

Overview of Activities for Next Week



In Lab Next Week:

1. Determine total mass of dry NaX product
2. Acquire IR spectrum of solid NaX product
3. Perform cobalt-exchange reaction:




1. Percent Yield of NaX Product

You prepared two solutions – sources of alumina and silicate – and mixed them together to form NaX:



Overall reaction:



How will you calculate the theoretical yield of NaX?

Find limiting reactants from preparation of Solutions 1 & 2 and for the overall reaction