DATE	EXP. NUMBER	EXPERIMENT		
8/28/2012	EAC-14	Preparation of DEET		
NAME			LAB PARTNER	WITNESS
Eim A. Chemist				

<u>Purpose:</u> To prepare N, N-diethyl-m-toluamide (DEET) and characterize it by IR and proton NMR

## Reaction:

1. 
$$CH_3 + SOCI_2 \xrightarrow{A} CH_3 + SO_2 + HOI$$

2.  $CH_3 + 2(CH_3CH_2)NH \xrightarrow{Emen} CH_3CH_3 + CH_3CH_3$ 

$$CH_3 + CH_3CH_2 \xrightarrow{C} NH_2CH_3 + CH_3CH_3 +$$

Physical Data and Safety:

Reagent	Structure/Formula	Mol. Wt.	Mp	Вр	Density	Safety
m-Toluic Acid	CHD2	136,15	108 -112°C	236 °C	1.05 g/mL	Causes eye, skin, respiratory tract irritation
Thionyl Chloride	SOCI <sub>2</sub>	118.97	-105°C	76 °C	1,63 g/mL	May be fatal if inhaled, Causes burns by all exposure routes, Water-reactive, Contact with water liberates toxic gas, Lachrymator (substance which increases the flow of tears).
Sodium Sulfate	Na <sub>2</sub> SO <sub>4</sub>			_	_	Possible irritant
Diethylamine	(CH3CH2)2NH	73,13	-50 °C	55 °C	0,7 g/mL	Harmful if absorbed through the skin, Causes eye and skin burns, Causes digestive and respiratory tract burns, Extremely flammable liquid and vapor, Vapor may cause flash fire,
Diethyl Ether	(CH <sub>8</sub> CH <sub>2</sub> ) <sub>2</sub> D	74,12	-116 °C	35 °C	0.7 g/mL	Extremely flammable liquid and vapor. May form explosive peroxides. Irritating to eyes, skin and respiratory tract. Vapors may cause drowsiness and diziness.
Sodium Hydroxide, 10% aq.	NaOH					Corrosive. Causes eye and skin burns. May cause respiratory tract burns.

DATE	EXP. NUMBER	EXPERIMENT		
8/28/2012	EAC-14	DEET		
NAME			LAB PARTNER	WITNESS
EA Chemist				

Reagent	Structure/Formula	Mol. Wt.	Mp	Вр	Density	Safety
Hydrochloric Acid, 10%	HCl					Causes eye, skin, and respiratory tract burns . Corrosive.
N, N-Diethyl-m- toluamide	CH3	191,27		111°C	0,996 g/mL	Irritant
Alamina	Al <sub>2</sub> D <sub>3</sub>					May cause respiratory tract irritation. May cause lung damage. May cause mechanical eye and skin irritation.
Hexane	CH14	86,18	-95 °C	62 <b>-</b> 69 °C	0,678 g/mL	Flammable liquid and vapor, Harmfal by inhalation, in contact with skin and if swallowed, Irritating to eyes and skin. Vapors may cause drowsiness and diziness,


When heating a reaction apparatus, be sure that it is open to the air so that pressure build up and subsequent rupture of the apparatus does not occur.

When heating liquids, make sure the liquid is stirred (or a boiling chip is added) to prevent "bumping".

When performing an extraction, make sure to vent the separatory funnel often to prevent pressure build-up.

## <u>Procedure</u>

1. Assemble the apparatus sketched below, be sure to vent the outlet into the hood! (NOTE: Apparatus must be dry!)



2. Add 0.272 g m-toluic acid to the reaction flash.

DATE	EXP. NUMBER	EXPERIMENT		
8/28/2012	EAC-14	DEET		
NAMÉ			LAB PARTNER	WITNESS
EA Chemist				

- 3. In a HOOD(!), add 0.30 mL of thionyl chloride.
- 4. Start circulation of water through the reflux condenser, then heat the reaction gently until the mixture begins to boil.
- 5. Reflux the mixture gently for 15 min, then cool to room temperature.
- 6. Using the syringe, add 4.0 mL anhydrous ether to the mixture and stir until solution is obtained.
- 7. Dissolve 0.66 mL of ice-cold diethylamine in 1.33 mL of anhydrous diethyl ether
- 8. Add solution from step #7 via syringe to the reaction flask dropwise over 10 15 min.
- 9. Stir reaction mixture for 10 min at room temperature.
- 10. Add 2 mL 10% NaOH (via syringe) and stir reaction mixture for 15 min at room temperature.
- 11. Transfer to separatory funnel and discard aqueous (bottom) layer.
- 12. Wash organic layer with 2 mL 10% NaOH. Discard aqueous (bottom) layer.
- 13. Wash organic layer with 2 mL 10% HCl. Discard aqueous (bottom) layer.
- 14. Dry the organic layer with sodium sulfate, then decant to a small conical vial.
- 15. Evaporate ether using a hot plate on a low setting.
- 16. Filter the crude product from step 15 through 1.0 g alumina in Pasteur pipet column, elute with hexane (5 7 mL) into a small tared beaker.
- 17. Evaporate the hexane on a hot plate over low heat.
- 18. Weigh product, obtain IR and proton NMR spectra.

DATE	EXP. NUMBER	EXPERIMENT		]
8/28/2012	EAC-14	DEET		
NAME	<b>.</b>		LAB PARTNER	WITNESS
EA Chemist				

## Data and Observations

<u>Step</u>	Data/Observations/Calculations
2.	Wt of m-toluic acid: 0.275 g. 0.275 g (1 mol/136.15 g) = 0.0020 mol
3,	Vol SOCl <sub>2</sub> = 0.30 mL, 0.30 mL(1.63 g/1 mL) = 0.489 g,
	0489 g (1 mol/118,97 g) = 0,0041 mol
5.	Refluxed 20 min, rxn took about 10 min to cool to room temperature.
6.	4.3 mL ether added, solution was a brown color
7.	Vol diethylamine = 0,70 mL, 0,70 mL(0,7 g/1 mL) = 0,49 g
	0,49 g (1 mol/73,13 g) = 0,0067 mol
	Vol ether used = 1.5 mL
8.	During the addition, a thick white cloud formed over the reaction mixture.
10.	Vol 10% NaOH = 2,4 mL
12, 13	Vol 10% NaOH = 2 mL, vol 10% HCl = 2 mL
14.	Enough sodium sulfate was added so that free-flowing granules were present, solution allowed to stand over sodium sulfate for about 15 min.
75,	Hot plate setting = 3, slow stream of air was directed over the top of the vial to assist evaporation. Crude product was a dark-brown liquid.
16,	Alamina removed the brown color from the product,
17.	Hot plate setting = 3, slow stream of air was directed over the top of the vial to assist evaporation.
	Product was a yellow liquid

DATE	EXP. NUMBER	EXPERIMENT		
8/28/2012	EAC-14	DEET		
NAME			LAB PARTNER	WITNESS
EA Chemist				

18. Beaker tare wt = 31,873 g Beaker + DEET = 32,213 g Wt DEET = 0,340 g (pale yellow oil).

Theoretical yield calculation:

Limiting reagent is m-toluic acid (MTA, 0.0020 mol)
0.0020 mol MTA(1 mol DEET/1 mol MTA)(191.27 g DEET/1 mol DEET) =
0.383 g DEET

Percentage yield = (0.340 g/0.383 g) x 100 = 88.8%

FTIR (film, NaCl plates): 2980, 2880 (sp³ C-H), 1633 (C=O), 1585 (aromatic C=C) cm⁻¹ IR corresponds to that of an authentic sample (see Spectral Database for Organic Compounds)

 $^{1}$ H-NMR (CDCl $_{3}$ ):  $\delta$  7.3 - 7.1 (m, 4H, Ar-H), 3.53 (q, 2H, ethyl CH $_{2}$ ), 3.24 (q, 2H, ethyl CH $_{2}$ ), 2.35 (s, 3H, Ar-CH $_{3}$ ), 1.23 (t, 3H, ethyl CH $_{3}$ ) 1.09 (t, 3H, ethyl CH $_{3}$ ) NMR corresponds closely to that of an authentic sample (see SDBS)

## Conclusions

DEET could be produced in 88.8% yield from m-tolaic acid, and the identity was confirmed by IR and proton NMR